MAQUET

Axius™ Coronary Shunt

Axius™ Coronaire Shunt

Dérivation Coronaire Axius™

Axius™-Koronarshunt

Shunt Coronarico Axius™

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

Figure 1.

Axius™ Coronary Shunt

1. Coil End Shunt Tip
2. Non-Coil End Shunt Tip
3. Threaded Radiopaque Tab with Tip Size
4. Lumen

INSTRUCTIONS FOR USE

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions.

DESCRIPTION:

The *Axius™* Coronary Shunt is a polymer tube with an embedded stainless steel coil with atraumatic tips at each end. There are openings in the tips to allow for perfusion of blood through the shunt and beyond the arteriotomy. The *Axius™* Coronary Shunt sizes are based on the outer diameter of the polymer tips. A radiopaque tab which identifies the shunt diameter is attached by a thread to the shaft of the shunt to facilitate positioning and removal from the arteriotomy. The *Axius™* Coronary Shunt is available in various sizes to fit various size coronary arteries.

INDICATIONS:

The *Axius™* Coronary Shunt is designed to help reduce blood in the operative field by temporary occlusion of the artery and to provide blood flow distal to the arteriotomy during the anastomosis. The *Axius™* Coronary Shunt is not an implant and is removed prior to completion of the anastomosis.

WARNINGS AND PRECAUTIONS:

- Use care during insertion and removal of the shunt to avoid tearing of the vessel wall and/or intimal dissection.
- During insertion of the shunt, avoid engaging shunt in septal or side branches.
- Discontinue use and remove if significant resistance is felt during insertion of the shunt.
- Care should be exercised upon removal of the shunt to prevent the shunt from getting entangled in sutures.
- Take care not to stitch the shunt while suturing anastomosis.
- Once shunt has been removed, do not reinsert as residual blood may remain in the lumen.

INSTRUCTIONS FOR USE:

1. Prepare coronary arteriotomy as for standard bypass operation.
2. Select the *Axius™* Coronary Shunt according to the internal diameter of the vessel.
3. Using forceps, gently insert one end of shunt into the arteriotomy far enough such that the other end can then be easily inserted into the arteriotomy. Discontinue insertion if significant resistance is felt during insertion and select a smaller size shunt.
4. Use threaded tab or forceps to center shunt in arteriotomy.
5. If the appropriate shunt size has been selected, there should be minimal blood seeping from the arteriotomy. If there is excessive seepage, remove the shunt and select a more appropriate size. Before inserting a larger shunt, verify that a septal branch is not causing bleeding at the arteriotomy site.
6. Prior to completion of the anastomosis, gently remove the shunt by pulling on the tab. Forceps can be used if required. If shunt cannot be removed with ease, verify that the shunt has not been sutured in the anastomosis. If this has occurred, remove the anastomosis suture from the shunt, as necessary, and resuture the vessel.

Note: Flow rate information is available upon request.

HOW SUPPLIED:

The Axius™ Coronary Shunt is supplied **STERILE** and **NON-PYROGENIC** in an unopened and undamaged package. It has been sterilized using gamma irradiation and is for single use only. **DO NOT RESTERILIZE. DO NOT REUSE.**

CAUTION: Federal (U.S.A.) law restricts this device to sale distribution and use by or on the order of a physician.


Figure 2.
Insertion of Coil End Shunt Tip

Using forceps, gently insert the coil end of the shunt into the arteriotomy.


Figure 3.
Insertion, Non-Coil End of Shunt

Coil end of shunt is inserted into the arteriotomy far enough that the non-coil end can be inserted. Shunt can be centered after both ends are within the arteriotomy.



Figure 4.
Shunt Removal

Pull the threaded tab attached to the shunt. Shunt will exit the arteriotomy asymmetrically.

WARRANTY:

MAQUET Cardiovascular LLC (MAQUET) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond MAQUET's control directly affect the instrument and the results obtained from its use. MAQUET's obligation under this warranty is limited to the repair or replacement of this instrument and MAQUET shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. MAQUET neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **MAQUET assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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