



<b>Packaging:</b> 1 Sterile Bag 04187 1 Label big 04148 0.05 Label big 04148 0.05 Label big PCTS 1030284 0.05 Paper between 04265 0.05 Carton 04120 0.05 IFU Polystan 1020969	Country <b>ARG</b> Institution <b>Cardiomedic</b> City	POLYSTAN TUBING SET <small>This document contains confidential and proprietary information. All design, reproduction, manufacturing, use and sale rights in such information are reserved.</small> Article-No. <b>1040939</b> Rev.-No. <b>01</b> Catalog-No. <b>802150</b> Drawn by <b>R.Brandl</b> Date <b>07-09-27</b> Verified Page <b>1 / 1</b>
<b>1</b> 00711 /RE "01405" PVC-70Sh 1/4x1/16x95 cm 00282 1/4x1/4 "01336" Silicone 1/4x1/16x55 cm 00282 1/4x1/4 "05000" PVC-75Sh 4.33x1.3x20 cm T-Conn. 05502 "05000" PVC-75Sh 4.33x1.3x20 cm 06303 05300 LLf/RE 00601	<b>2</b> 00501 "05001" PVC-75Sh 3.5x1.0x 100 cm T-Conn. 05502 "05001" PVC-75Sh 3.5x1.0x 100 cm 00306 LLf 00833 /WH 00402 LLf 00833 /WH no cap "05001" PVC-75Sh 3.5x1.0x 5 cm 00171 Trans-Protector GVS TP 10 "05000" PVC-75Sh 4.33x1.3x170 cm 00307 LLm 00601	1 Sterile Bag 04187 1 Sterile Bag 04136

